

optional ATC Temp. probe for pH meter

TEMPERATURE PROBE

Model : TP-06, TP-07

ISO-9001, CE, IEC1010



FEATURE :

Model : TP-06

- * The optional Temp. probe (ATC probe) for pH meter PH-206.
- * Transistor sensor.
- * Range : 0 to 60 °C.
- * Accuracy : ± 0.8 °C/1.5 °F.

Model : TP-07

- * The optional Temp. probe (ATC probe) for pH meter PH-207, PH-207HA, PH-208, PH-221, YK-2005WA, WA-2017SD, PH-221E, BPH-231.
- * Thermistor sensor.
- * Range : 0 to 60 °C.
- * Accuracy : ± 0.8 °C/1.5 °F.

Lutron

LUTRON ELECTRONIC



The Art of Measurement